



Link **SymphoKnee**
Built on **TRUST.**

LinkSymphoKnee

Built on **TRUST.**



System Description:

- 02 Introduction
- 03 Implants Combination Overview
- 04 System Overview
 - *LinkSymphoKnee* CR
 - *LinkSymphoKnee* UC
 - *LinkSymphoKnee* PS / PS+
 - *LinkSymphoKnee* CCK

Design Philosophy

- 06 Design Philosophy
- 08 Sizes
- 09 Stems

Features and benefits:

- 10 Femoral Components CR/PS
- 12 Femoral Components CCK
- 14 Articulating Surfaces CR/UC
- 16 Articulating Surfaces PS/PS+/CCK
- 18 Tibial Components Monoblock/Modular
- 20 Augments
- 22 Stems
- 26 LINK PorEx (TiNbN)
- 27 E-Dur
- 28 Surgical Technique
- 28 Instruments

Sizes/Dimensions

- 30 Femoral Components CR/PS/CCK
- 31 Tibial Components Monoblock/Modular
- 31 Patella Components
- 32 Articulating Surfaces CR/UC/PS/PS+/CCK
- 34 Cylindrical Press-Fit Stems
- 35 Cylindrical Cementless Stems
- 36 Conical Cementless Stems
- 37 Conical Cemented Stems
- 37 Centralizers
- 38 Femoral Augments
- 39 Tibial Augments

Literature

Important Information

LinkSymphoKnee

Introduction

Rooted in trust and decades of expertise, the *LinkSymphoKnee* sets a new standard in bicondylar knee joint prostheses. Its design integrates knowledge from proven implant systems, advanced kinematic concepts, and cutting-edge material and coating technologies, ensuring optimal performance and adaptability.

The *LinkSymphoKnee* is a versatile, complete knee system that supports a wide range of treatment scenarios, from primary procedures to complex revision cases.

Key Features and Configurations

- Flexibility in design:**

The system offers configurations that allow for posterior cruciate ligament (PCL) retention or substitution, as well as a condylar constraint knee (CCK) option. These configurations can be enhanced with stems and augments for additional support in more demanding cases.

- Available configurations: **CR** (Posterior Cruciate ligament-**R**etaining), **UC** (**U**ltra **C**ongruent), **PS** (Posterior Cruciate ligament-**S**acrificing), **PS+** (Posterior Cruciate ligament-**S**acrificing with varus/valgus stability), **CCK** (Condylar **C**onstrained **K**nee)

- Tibial component options:**

The tibial component is available in two designs:

- Monoblock: A single-piece option for simplicity and durability, including an All-Poly PS tibial variant.
- Modular: Compatible with all *LinkSymphoKnee* configurations, offering enhanced adaptability with stems and augments.

- Material excellence:**

Components are crafted using high-quality materials for durability and biocompatibility:


- EndoDur (CoCrMo) or LINK PorEx (CoCrMo + TiNbn HighGliss) for superior wear resistance.
- UHMWPE or E-Dur (HXLPE + Vitamin E) for improved polyethylene performance.

- Fixation versatility:**

The *LinkSymphoKnee* supports both cemented and uncemented fixation methods, offering surgeons flexibility to meet the needs of diverse patient profiles.


Designed to adapt to the complexities of modern knee arthroplasty, the *LinkSymphoKnee* combines innovation, reliability, and precision to deliver excellent patient outcomes. Whether for primary cases or complex revisions, this system is engineered to restore mobility and enhance quality of life.

LinkSymphoKnee



LinkSymphoKnee

Implants Combination Overview



* Uncemented components are not available for FDA markets
 ** LinkSymphoKnee UC (Ultra Congruent) released for FDA markets only

LinkSymphoKnee

System Overview



LinkSymphoKnee CR

(Posterior Cruciate ligament-Retaining)

- Requires intact PCL
- Requires stable collateral ligaments
- The Femoral component same as UC

LinkSymphoKnee UC

(Ultra Congruent)

- PCL can be sacrificed, whilst femoral bone stock is preserved.
 - Can optionally be used as a CR solution
- Requires stable collateral ligaments.
- Femoral component same as CR



LinkSymphoKnee PS

(Posterior Cruciate ligament-Sacrificing)

- PCL is sacrificed
- Requires stable collateral ligaments
- Femoral component same as PS+

LinkSymphoKnee PS+

(Posterior Cruciate ligament-Sacrificing with varus/valgus stability)

- PCL is sacrificed
- Requires collateral ligaments
- Femoral component same as PS.



LinkSymphoKnee CCK


(Condylar Constrained Knee)

- PCL is sacrificed
- Requires collateral ligaments
 - Provides implant stability where soft tissue balance is not satisfying
- Stems can be combined with modular tibia and CCK femur
- Femoral component also compatible with the PS and PS+

LinkSymphoKnee


Design Philosophy

The *LinkSymphoKnee* is a bone-preserving posterior-referenced system with easy flexion gap balancing.



Narrow femoral component design


This ensures better position and alignment of the implant. For wider native femora the LSK range comprises “+” sizes which increase the dimension in M-L of the implant.



- **Narrow design:** Prevents overhang, ensuring optimal fit for better patient outcomes.
- **Soft tissue protection:** Minimizes the risk of soft tissue impingement, reducing complications.
- **Improved longevity and comfort:** Enhances implant survival rates and reduces patient pain for better long-term results.¹
- **Versatile alignment:** Accommodates a broader range of external rotation in the native knee, providing greater flexibility in treatment options.
- **Simplified sizing:** Eliminates the need to match native external rotation to select the largest femoral component, streamlining the surgical process.²

Condyle design

- The *LinkSymphoKnee* features a condyle design with a **decreasing multi-radius**.
- **Key benefits of this design:**
 - **Optimized stability:** Effectively manages the shifting forces in the knee during movement.
 - **Smooth motion transition:** Ensures seamless movement from flexion to extension.
 - **Reduced midflexion instability:** Minimizes the instability often experienced in mid-flexion, enhancing patient confidence and joint performance.
- This innovative design supports better knee functionality and improved outcomes during dynamic motion.





Patella-friendly Design


Anterior Knee Pain (AKP) and the *LinkSymphoKnee* solution

- **AKP in TKA:**
 - A leading cause of early revision surgery, often caused by offset errors, oversizing, rotational errors, or patellar maltracking.^{4,5} In TKA the patella is often medialized by 2.5 mm compared to the native knee.⁶

- **Patella-friendly *LinkSymphoKnee* design:**
 - **Prevents maltracking:** Wide, shallow and lateralized trochlea groove optimizes patellar alignment.
 - Accommodates a variety of Q-angles.
 - The S-curve on the patellofemoral track ensures a smooth transition during flexion/extension.
 - **Minimizes offset errors:**
 - Thin patella shield prevents overstuffing.
 - 4.5° angled shield reduces notching and improves cement layer compression.
 - **Narrow femoral design:** Matches A/P dimensions without M/L overhang.
 - **Surgical versatility:** The option to trial the final patella shield position by using a drill-pin before the bone is cut ensures a precise fit of the implant.
 - This design addresses AKP risks, enhancing fit, alignment, and long-term outcomes.



Trochlea groove widens in extension



lateral medial

Sizes




LinkSymphoKnee – Femoral Components

- 14 femoral sizes (CR, PS, CCK)
- Materials: EndoDur (CoCrMo) and LINK PorEx (TiNbN)



LinkSymphoKnee – Tibial Components

- 10 tibial sizes, monoblock as well as modular
- Materials: EndoDur (CoCrMo) and LINK PorEx (TiNbN)



LinkSymphoKnee – Articulating Surfaces

- 6 insert heights for CR, UC and PS/PS+ (10–18 mm) for each tibia size
- 8 insert heights for CCK (10-24 mm) for each tibia size
- Materials: UHMWPE or **e-dur**



LinkSymphoKnee – All-Poly PS Tibial Components

- 3 insert heights (10, 12 and 14 mm) for All-Poly PS Tibial Components for each tibia size
- Material: UHMWPE

Stems

Conical Cemented Stems

- 10 length options (50 – 280 mm length)
- 4 diameters (10 – 15 mm)
- 3 offsets (0 mm, 3 mm, 6 mm)
- Material: EndoDur (CoCrMo) and LINK PorEx (TiNbN)



Cylindrical Press-Fit Stems

- 6 length options (80 – 240 mm)
- 11 diameters (10 – 20 mm)
- 3 offsets (0 mm, 3 mm, 6 mm)
- Material: Tilastan-S (Ti6Al4V)



Cylindrical Cementless Stems

- 5 length options (80 – 240 mm)
- 11 diameters (10 – 20 mm)
- 3 offsets (0 mm, 3 mm, 6 mm)
- Material: Tilastan-S (Ti6Al4V)



Conical Cementless Stems

- 3 length options (80 – 240 mm)
- 13 diameters (10 – 20 mm)
- No offset options
- Material: Tilastan-S (Ti6Al4V)



Femoral and Tibial Augments

- 10 sizes
- 3 heights (5 – 15 mm)
- Material: Tilastan-S (Ti6Al4V)



LinkSymphoKnee

Femoral Components CR/PS 



Anatomical notch geometry:
Cruciate-friendly, patella-friendly




Range of sizes

- 14 femoral sizes
- The femoral component growth evenly in 3 mm A-P



Mechanism:
The pegs on the CR provide initial stability and orientation of the femur

Groove depth
Reduced risk of PF ligament strain⁶



Equal distal and posterior thickness of 9 mm
Facilitates extension and flexion gap consistency




Material:
Available in EndoDur (CoCrMo) or LINK PorEx (TiNbN)



Fixation type:
The *LinkSymphoKnee* is available in cemented and uncemented fixation in combination with LINK PorEx only





Mechanism:

- Improved jump height
- Less paradoxical motion
- Improved kinematics
- Early post-cam engagement

Concave cam with increased contact surface

Bone-preserving
Angled and open PS box

Narrow anterior flange
Forgiving to soft tissue

Wide proximal trochlear groove

Medialized trochlear groove

Graduate coronal curvature
Larger sizes feature a wider gap between contact patches


LinkSymphoKnee

Femoral Components CCK

Modular Stems:

Increases stability in poor bone conditions⁷

- Cemented (EnduDur-S, LINK PorEx):
- Press-fit (Tilastan-S)
- Cementless, conical (Tilastan-S)
- Cementless, cylindrical (Tilastan-S)



Fixationtype:

Cemented or cementless fixation with stems for additional stability




Range of sizes

- 14 femoral sizes
- The femoral component growth evenly in 3 mm A-P



Material:

Available in LINK PorEx (TiN**b**N) for each *LinkSymphoKnee* configuration





Reduced risk of overstuffing

Minimized patellar flange thickness

Anatomic stem ML and AP location
 Reduced need for offset stems due to anatomical position of the stem taper connection


Reduces mid-flexion instability
 Gradually decreasing flexion radius



Mechanism:
 The stem / Taper Cap helps to align the Femoral Component, provides initial stability and neutralizes rotational forces.⁷



Reduced risk of notching
 4.5° patellar flange angle



Equal distal and posterior thickness of 9 mm
 Facilitates extension and flexion gap consistency


Stem connection and stem offset start close to joint line

LinkSymphoKnee

Articulating Surfaces CR/UC



Material:
Available in UHMWPE or E-Dur (Vitamin E infused highly cross-linked UHMWPE)



2up / 2down compatibility
Flat mediolateral edge at posterior facet

Soft tissue sparing
High-flexion recess for patellar ligament

LinkSymphoKnee – Articulating Surface CR

Posterior lip:
0.7-1.0 mm
size-dependent

Anterior lip:
2.6-5.4 mm
size-dependent


Optimized for 3° bone resection



Triply-secured inserts:
AP and ML self-locking press-fit
with antirotational guiding bar





Fixation type:
Dovetail locking mechanism



• **Dovetail locking mechanism**


The dovetail locking mechanism works via press-fit and does not require any additional screw. This locking mechanism is used for CR, UC, PS, PS+ and CCK Articulating Surfaces.

Soft tissue sparing
Wide intercondylar notch



2up / 2down compatibility
Flat mediolateral edge at posterior facet

LinkSymphoKnee – Articulating Surface UC




LinkSymphoKnee

Articulating Surfaces PS/PS+/CCK



Material:
Available in UHMWPE or E-Dur (Vitamin E infused highly cross-linked UHMWPE)




PS jump height
17 mm

PS+ jump height
17 mm

PS Constraint
Free to rotate




Recess for high flexion capability
Patella tendon clearance



Fixation type:
Dovetail locking mechanism

PS+ Constraint
Constrained 5° IE and 3° V/V

Posterior Chamfer
reduces impingement



Hyperextension

PS: 12°
PS+: 10°
CCK: 10°

The differences between the posterior-stabilized inserts


CCK jump height
17 mm

Recess on CCK
post for high
flexion capability

- PS
- PS+
- CCK

Varus/Valgus (v/v) constrained			
	PS	PS+	CCK
IE Rotation	not limited	±5°	±2.5°
Varus/Valgus	±4°	±3°	±1°

16.5 mm
17.3 mm
17.4 mm



CCK Constraint
Constrained 2.5° IE and 1° VV



Introperative flexibility



Small inventory many options (compatibility)

- **Femur insert matching**
To allow the best fit for the patient, the femur is compatible with two sizes up and down on the tibia.

Femoral Size 5



Insert Size 3-4



Insert Size 5-6



Insert Size 7-8



Tibia Size

3 4



Tibia Size

5 6




Tibia Size

7 8

Bonus

LinkSymphoKnee

Tibial Components




Material:
Available in EndoDur (CoCrMo) or LINK PorEx (TiNbn)




Fixationtype:
The *LinkSymphoKnee* is available in cemented and uncemented fixation in combination with LINK PorEx only





High flexibility & low storage need


The dovetail locking mechanism is the same on all inserts. This allows for a small inventory of tibial components despite a variety of treatment options.



LinkSymphoKnee

Augments


The *LinkSymphoKnee* offers a range of femoral and tibial augments designed to address bone deficiencies. Femoral augments compensate for distal and posterior femoral defects and pair with CCK femoral components, while tibial augments address proximal tibial defects and pair with modular tibial components. All screws are preassembled, ensuring a swift workflow. The tibial augments are conical and reduce the distal resection interface by up to **3** sizes tibially.



Material:
Tilastan-S (Ti6Al4V)



LinkSymphoKnee – Tibial Augments






Wide range of sizes:
5 - 15 mm in height



Fixation:
1 pre-assembled screw
tightened with 3 Nm



Efficient surgical technique:
No additional instrument set needed

LinkSymphoKnee – Femoral Augments

MEDIAL AUGMENTS

LATERAL AUGMENTS

Posterior Options

Medial Augment Size	Posterior Option Size	Lateral Augment Size
15 mm	5 mm	15 mm
10 mm	5 mm	10 mm
5 mm	5 mm	5 mm

LinkSymphoKnee

Stems

What is the function of stems in the TKA?

- Additional implant fixation in the diaphysis
- Load transfer from compromised metaphyseal bone structure to intact diaphysis bone structure.

“Although the use of the stem remains controversial in constrained TKA, the stem should be used routinely if there is inadequate bony surface.”¹⁰

Rotation

Possibility of 360° of rotation for offset stems




Markings

For orientation of offset stems



Easy to revise

The taper adapter design disengages the stem from the modular tibia or CCK femur should a need for revision arise.



Wide Variety

The *LinkSymphoKnee* offers an extensive stem portfolio, from 0 mm offset straight cemented up to 6 mm offset press-fit stems.



Material

- Cemented: EndoDur (CoCrMo) and LINK PorEx (TiNbN)
- Cementless: Tilastan-S






High flexibility & low storage need

- **Stem femur/tibia matching**

To minimize storage space while ensuring best-fit solutions, the stem portfolio has only a female coupling and is compatible with the femur as well as the tibia component of the *LinkSymphoKnee*.



Offset Magnitude

0 mm, 3 mm and 6 mm

Why do we need offset

The significant variability of the bone canal relative to the femoral or tibial component necessitates adaptation of longer stems with an offset.

The femoral canal is often located anterior to the center of the femur component.

The tibia canal is often located anterior and medial to the plateau center.¹¹



Wide range of sizes

- Cemented: 10 length (50 - 280 mm) and 4 diameters (10, 11, 13, 15 mm)
- Cementless: up to 6 length (80 - 240 mm) and 12 diameters (10 - 20 in 1-mm steps and 22 mm)



Centralizers

- The cemented stem must be assembled with a centralizer to find the ideal seating in the canal.

Allow a central position of the cemented stem in the medullary canal.


This prevents contact between the metal stem and the corticalis and therefore stress peaks in the bone in case of bending load.⁶

LinkSymphoKnee

Stems

Cylindrical Press-Fit Stems

- The **slotted stem design** with **rib structure** at the tip enhances **rotational stability**, prevents stress shielding, and reduces stem stiffness. This helps lessen **stem tip pain** in patients.¹²
- The **minimal roughness** of the stem is intended to provide a **firm grip** without full bone integration.
- Reduces the risk of **puncturing** or **fracturing** the bone during implantation of a **cementless stem**.¹³
- The combination of a **slotted press-fit design** and the **minimal roughness** of the stem surface allow **bone-preserving revision** if needed.



0.6-mm built-in press-fit



0.6-mm built-in press-fit



Cylindrical Cementless Stems

- The **cylindrical press-fit stems** provide **diaphyseal anchoring** with osteointegration.
- This design distributes the **load** from compromised metaphyseal bone to the intact **diaphysis**.
- **Ribs** forming a **star geometry** around the stem ensure **rotational stability** while distributing loads evenly and minimizing **peak forces**.
- The combination of ribs and the **PoroLink surface** promotes **bone ongrowth**, anchoring the stem firmly to handle **rotational, push, and pull forces**.¹⁴
- Highly effective in cases involving **periprosthetic fractures**.¹⁵



Efficient surgical technique


One small tray with only 6 additional instruments




Fixation Type

Cemented and cementless

1-mm
cement mantle



2.5-mm
cement mantle




Conical Cemented Stems

- **Cemented stems** are made from **CoCrMo** and are also available with a **LINK PorEx coating**.
- Offering **freedom of placement** with the **short stem version**.
- The combination of **round** and **flat geometries** provides **rotational stability** while minimizing **peak stress** at the **stem/cement interface**.
- Manufactured from the same material as the implant.
- The **LINK PorEx coating** reduces **nickel ion release** and enhances **tribological properties**.¹⁶
- Used with the **Endo-Model Knee System** since 1992 with excellent results.

Conical Cementless Stems

- The **conical press-fit stem** features ribs that increase the **surface area** for anchoring and distribute the **load** across the **metaphysis** and **diaphysis**.
- **Optimized Press-Fit:** As the stem diameter increases the press-fit decreases, allowing for **precise positioning**.
- The **conical design** allows for **self-centering preparation**, minimizing the risk of stem tilt.
- Used with the **Endo-Model Knee System** since 1992, delivering consistent success.



LinkSymphoKnee – LINK PorEx (TiNbN)

Surface coatings:

The *LinkSymphoKnee* family with PorEx has a ceramic-like surface.¹⁷




Material:

Available in LINK PorEx (TiNbN) for each *LinkSymphoKnee* configuration



Wettability:

LINK PorEx has high wettability, supporting lubrication in the joint.¹⁷



Hardness:

LINK PorEx has outstanding hardness.¹⁷



Metal ion release:

Significant reduction of metal ion release.¹⁷



ION



Surface coatings:

Protect the sliding surface from wear and scratches.¹⁷

ε-dur[®] Description

As early as the 1960s, studies confirmed that UHMWPE is an outstanding material for sliding bearings in joint replacements.^{18,19,20}

The key properties of UHMWPE are:


- Biocompatibility
- High abrasion resistance^{19,20}
- Impact resistance
- Fatigue strength and crack resistance

It therefore conforms to national and international standards for implant materials.^{18,20}


Consequently, UHMWPE became the “gold standard” for tibial articular surfaces in knee joint reconstructions.^{19,22} This status is confirmed by

outstanding long-term results of published studies and registry data.^{23,24,25} Now this “gold standard” for orthopedic implants has been improved by further increasing its mechanical properties and durability.²⁶ High crosslinking produced a marked improvement in abrasion resistance.²⁷⁻³² In addition, effective oxidation protection is achieved by enriching the material with vitamin E.^{27,28}

In the production of ε-dur[®] Vit-E plateaus, vitamin E is used as an antioxidant in order to protect the material by neutralizing the free radicals created by highly cross-linking. The product’s mechanical properties and biocompatibility are preserved.^{19,27,33-35}



Properties




Wear Resistance³⁶

The process of intensive cross-linking improves the abrasion resistance of UHMWPE, while enrichment with vitamin E counteracts the aging process.


LinkSymphoKnee – Surgical Technique

The *LinkSymphoKnee* instruments follow the surgeon's preferences allowing for extra-medullary or intra-medullary alignment of the tibia as well as the option for femur first, tibia first, or distal cut first.



LinkSymphoKnee – Instruments

The *LinkSymphoKnee* instruments are a streamlined modular system designed to provide accurate and reproducible results, with a minimized inventory. In a primary CR or PS case only three instrument trays are needed. Each further option, like a macro/micro or wide size, requires just one additional tray.



For revision surgery with the *LinkSymphoKnee* CCK, a minimum of only four instrument trays are needed. The intuitive and elegant instruments are inspired by the SPAR-K instruments from the Gemini SL. With the golden color it is easy to recognize moving parts and where to adjust the instruments.

Distal Cut Resection Guide

Accurately determines bone resection level even for severe defects with the option to adjust the varus/valgus.



Femoral Sizing Guide

The femur size is determined using the posterior condyles as a reference. The femur rotation can be adjusted with the golden knob.

A condylar defect can be addressed with magnetic shims.



Offset Selector Guide

The CCK stem offset selector is a handy instrument to determine the correct offset if a stem is being used. This instrument can be used for left and right knees and has built in 6° valgus. The golden knob can be turned 360° while representing 0 mm, 3 mm, or 6 mm offset.

Feature: The determined offset and rotation is “saved” on the instrument for later stem and implant coupling.




Rotation:
Possibility of
360° of rotation




Offset Magnitude:
0 mm, 3 mm and 6 mm

LinkSymphoKnee Femoral Components – CR/PS




LinkSymphoKnee Femoral Components – CCK



Size	A	B	C/D	E	F	G	H	I
	AP mm	ML mm	Post & Distal Thickn. mm	Posterior Condyle Length mm	Condyle Width mm	Inner Box Width mm	Outer Box Width mm	Stem Location A-P mm
	CR/PS/CCK	CR/PS/CCK	CR/PS/CCK	CR/PS/CCK	CR	PS/CCK	PS/CCK	CCK
0	47	53.3	9	28.5	14.5	14.5	18.5	12
1	50	56	9	29.8	15.2	14.5	18.5	12.8
2	53	58.5	9	30.9	15.9	14.5	18.5	13.5
3	56	61	9	32.5	16.5	17.5	21.5	12.4
3+	56	63.5	9	32.5	16.5	17.5	21.5	12.4
4	59	63.5	9	34	17.2	17.5	21.5	13.6
4+	59	66	9	34	17.2	17.5	21.5	13.6
5	62	66	9	35.8	17.9	17.5	21.5	14.5
5+	62	69	9	35.8	17.9	17.5	21.5	14.5
6	65	69	9	37.3	18.7	17.5	21.5	14.8
7	68	72	9	38.9	19.5	17.5	21.5	15.5
8	71	75	9	40.1	20.3	17.5	21.5	16
9	74	78	9	41	21.1	17.5	21.5	17
10	77	81	9	42	22	17.5	21.5	17.5


LinkSymphoKnee Tibial Components




Size	A	B	C	D	E	F
	AP mm	ML mm	Plate Thickness mm	Stem Length mm	Stem Length mm	Slope
1	59	37.5	4	38.5	43.8	3°
2	62.5	40	4	38.5	43.8	3°
3	66	42.5	4	38.5	43.8	3°
4	69.5	45	4	38.5	43.8	3°
5	73	47.5	4	38.5	43.8	3°
6	76.5	50	4	38.5	43.8	3°
7	80	52.5	4	38.5	43.8	3°
8	83.5	55	4	38.5	43.8	3°
9	87	57.5	4	38.5	43.8	3°
10	90.5	60	4	38.5	43.8	3°

LinkSymphoKnee Patella Components

A	B	F
Ø mm	Height mm	Peg Height mm
25	6	5
28	6	5
31	7	5
34	8	5
37	9	5
40	10	5



LinkSymphoKnee Articulating Surfaces – CR, UC, PS, PS+, CCK



Size	A	B	C		D		E			F		G	H
	AP mm	ML mm	Post. Lip Height mm		Ant. Lip Height mm		Post Width mm			Post Height mm		Post. Lip Height mm	Ant. Lip Height mm
	All	All	CR	UC	CR	UC	PS	PS+	CCK	PS/PS+	CCK	All	All
1-2	59	37.5	0.8	3.5	3.1	5.5	13.5	14.3	14.35	18.4	21	2.1	2.1
1up	59	37.5	–	–	–	–	16.5	17.3	17.35	19.4	22	1.6	1.6
3-4	66	42.5	0.9	3.5	4.6	7.5	16.5	17.3	17.35	19.4	22	2.3	2.9
3-4down	66	42.5	–	–	–	–	13.5	14.3	14.35	18.4	21	2.9	3.9
5-6	73	47.5	1	3.5	5.6	9	16.5	17.3	17.35	20.4	23	3.1	3.6
7-8	80	52.5	1.1	3.5	6.6	9.5	16.5	17.3	17.35	21.4	24	3.5	4.6
9-10	87	57.5	1.2	3.5	7.4	11	16.5	17.3	17.35	22.4	25	4	5.1



Colour coding:
for streamlined workflow

The following table shows the possible *LinkSymphoKnee* – CR size combinations:

		CR Femoral Component											
		0	1	2	3/3+	4/4+	5/5+	6	7	8	9	10	
Tibial Component	1	CR Articulating Surface 1-2					x	x	x	x	x	x	x
	2	CR Articulating Surface 1-2					x	x	x	x	x	x	x
	3	x	CR Articulating Surface 3-4						x	x	x	x	
	4	x	CR Articulating Surface 3-4						x	x	x	x	
	5	x	x	x	CR Articulating Surface 5-6						x	x	
	6	x	x	x	CR Articulating Surface 5-6						x	x	
	7	x	x	x	x	x	CR Articulating Surface 7-8						
	8	x	x	x	x	x	CR Articulating Surface 7-8						
	9	x	x	x	x	x	x	x	CR Articulating Surface 9-10				
	10	x	x	x	x	x	x	x	CR Articulating Surface 9-10				

The following table shows the possible *LinkSymphoKnee* – UC size combinations:

		UC Femoral Component											
		0	1	2	3/3+	4/4+	5/5+	6	7	8	9	10	
Tibial Component	1	UC Articulating Surface 1-2					x	x	x	x	x	x	x
	2	UC Articulating Surface 1-2					x	x	x	x	x	x	x
	3	x	UC Articulating Surface 3-4						x	x	x	x	
	4	x	UC Articulating Surface 3-4						x	x	x	x	
	5	x	x	x	UC Articulating Surface 5-6						x	x	
	6	x	x	x	UC Articulating Surface 5-6						x	x	
	7	x	x	x	x	x	UC Articulating Surface 7-8						
	8	x	x	x	x	x	UC Articulating Surface 7-8						
	9	x	x	x	x	x	x	x	UC Articulating Surface 9-10				
	10	x	x	x	x	x	x	x	UC Articulating Surface 9-10				

The following table shows the possible *LinkSymphoKnee* – PS/PS+ size combinations:


		PS Femoral Component											
		0	1	2	3/3+	4/4+	5/5+	6	7	8	9	10	
Tibial Component	1	PS/PS+ Articulating Surface 1-2			PS/PS+ Articulating Surface 1-2up		x	x	x	x	x	x	x
	2	PS/PS+ Articulating Surface 1-2			PS/PS+ Articulating Surface 1-2up		x	x	x	x	x	x	x
	3	x	PS/PS+ Articulating Surface 3-4down		PS/PS+ Articulating Surface 3-4				x	x	x	x	
	4	x	PS/PS+ Articulating Surface 3-4down		PS/PS+ Articulating Surface 3-4				x	x	x	x	
	5	x	x	x	PS/PS+ Articulating Surface 5-6						x	x	
	6	x	x	x	PS/PS+ Articulating Surface 5-6						x	x	
	7	x	x	x	x	x	PS/PS+ Articulating Surface 7-8						
	8	x	x	x	x	x	PS/PS+ Articulating Surface 7-8						
	9	x	x	x	x	x	x	x	PS/PS+ Articulating Surface 9-10				
	10	x	x	x	x	x	x	x	PS/PS+ Articulating Surface 9-10				

The following table shows the possible *LinkSymphoKnee* – CCK size combinations:

		CCK Femoral Component											
		0	1	2	3/3+	4/4+	5/5+	6	7	8	9	10	
Tibial Component	1	CCK Articulating Surface 1-2			CCK Articulating Surface 1-2up		x	x	x	x	x	x	x
	2	CCK Articulating Surface 1-2			CCK Articulating Surface 1-2up		x	x	x	x	x	x	x
	3	x	CCK Articulating Surface 3-4down		CCK Articulating Surface 3-4				x	x	x	x	
	4	x	CCK Articulating Surface 3-4down		CCK Articulating Surface 3-4				x	x	x	x	
	5	x	x	x	CCK Articulating Surface 5-6						x	x	
	6	x	x	x	CCK Articulating Surface 5-6						x	x	
	7	x	x	x	x	x	CCK Articulating Surface 7-8						
	8	x	x	x	x	x	CCK Articulating Surface 7-8						
	9	x	x	x	x	x	x	x	CCK Articulating Surface 9-10				
	10	x	x	x	x	x	x	x	CCK Articulating Surface 9-10				


Symbol Description: x = Combination not permitted

Cylindrical Press-Fit Stems




Offset magnitude mm	D1 mm	L mm	LF mm	LT mm
0/3/6	10	80	115	98
0/3/6	11	80	115	98
0/3/6	12	80	115	98
0/3/6	13	80	115	98
0/3/6	14	80	115	98
0/3/6	15	80	115	98
0/3/6	16	80	115	98
0/3/6	17	80	115	98
0/3/6	18	80	115	98
0/3/6	10	120	155	138
0/3/6	11	120	155	138
0/3/6	12	120	155	138
0/3/6	13	120	155	138
0/3/6	14	120	155	138
0/3/6	15	120	155	138
0/3/6	16	120	155	138
0/3/6	17	120	155	138
0/3/6	18	120	155	138
0/3/6	19	120	155	138
0/3/6	20	120	155	138
0/3/6	11	160	195	178
0/3/6	12	160	195	178
0/3/6	13	160	195	178
0/3/6	14	160	195	178
0/3/6	15	160	195	178
0/3/6	16	160	195	178
0/3/6	17	160	195	178
0/3/6	18	160	195	178
0/3/6	19	160	195	178
0/3/6	20	160	195	178
0/3/6	12	200	235	218
0/3/6	14	200	235	218
0/3/6	16	200	235	218
0/3/6	18	200	235	218
0/3/6	20	200	235	218
0/3/6	22	200	235	218
0/3/6	12	220	255	238
0/3/6	14	220	255	238
0/3/6	16	220	255	238
0/3/6	18	220	255	238
0/3/6	20	220	255	238
0/3/6	22	220	255	238
0/3/6	12	240	275	258
0/3/6	14	240	275	258
0/3/6	16	240	275	258
0/3/6	18	240	275	258
0/3/6	20	240	275	258
0/3/6	22	240	275	258

Cylindrical Cementless Stems




Offset magnitude mm	D1 mm	D2 mm	L mm	LF mm	LT mm
0/3/6	16	10	80	115	98
0/3/6	16	11	80	115	98
0/3/6	16	12	80	115	98
0/3/6	16	13	80	115	98
0/3/6	16	14	80	115	98
0/3/6	16	15	80	115	98
0/3/6	16	16	80	115	98
0/3/6	17	17	80	115	98
0/3/6	18	18	80	115	98
0/3/6	16	10	120	155	138
0/3/6	16	11	120	155	138
0/3/6	16	12	120	155	138
0/3/6	16	13	120	155	138
0/3/6	16	14	120	155	138
0/3/6	16	15	120	155	138
0/3/6	16	16	120	155	138
0/3/6	17	17	120	155	138
0/3/6	18	18	120	155	138
0/3/6	19	19	120	155	138
0/3/6	20	20	120	155	138
0/3/6	16	11	160	195	178
0/3/6	16	12	160	195	178
0/3/6	16	13	160	195	178
0/3/6	16	14	160	195	178
0/3/6	16	15	160	195	178
0/3/6	16	16	160	195	178
0/3/6	17	17	160	195	178
0/3/6	18	18	160	195	178
0/3/6	19	19	160	195	178
0/3/6	20	20	160	195	178
0/3/6	16	12	200	235	218
0/3/6	16	14	200	235	218
0/3/6	16	16	200	235	218
0/3/6	18	18	200	235	218
0/3/6	16	12	240	275	258
0/3/6	16	14	240	275	258
0/3/6	16	16	240	275	258
0/3/6	18	18	240	275	258

Conical Cementless Stems



D1 mm	D2 mm	D3 mm	L mm	LF mm	LT mm
16	12	9	128	163	146
16	13	10	128	163	146
16	14	11	128	163	146
16	15	12	128	163	146
16	16	13	128	163	146
17	17	14	128	163	146
18	18	15	128	163	146
19	19	16	128	163	146
20	20	17	128	163	146
21	21	18	128	163	146
22	22	19	128	163	146
23	23	20	128	163	146
24	24	21	128	163	146
16	12	9	158	193	176
16	13	10	158	193	176
16	14	11	158	193	176
16	15	12	158	193	176
16	16	13	158	193	176
17	17	14	158	193	176
18	18	15	158	193	176
19	19	16	158	193	176
20	20	17	158	193	176
21	21	18	158	193	176
22	22	19	158	193	176
23	23	20	158	193	176
24	24	21	158	193	176
16	12	9	188	223	206
16	13	10	188	223	206
16	14	11	188	223	206
18	15	12	188	223	206
16	16	13	188	223	206
17	17	14	188	223	206
18	18	15	188	223	206
19	19	16	188	223	206
20	20	17	188	223	206
21	21	18	188	223	206
22	22	19	188	223	206
23	23	20	188	223	206
24	24	21	188	223	206

Conical Cemented Stems




Offset magnitude mm	D1 mm	D2 mm	D3 mm	L mm	LF mm	LT mm
0	16	11	8	50	85	68
0/3/6	16	10	7	80	115	98
0/3/6	16	11	8	80	115	98
0/3/6	16	13	10	80	115	98
0/3/6	16	15	12	80	115	98
0	16	10	7	95	130	113
0	16	11	8	95	130	113
0	16	13	10	95	130	113
0	16	15	12	95	130	113
0/3/6	16	10	7	120	155	138
0/3/6	16	11	8	120	155	138
0/3/6	16	13	10	120	155	138
0/3/6	16	15	12	120	155	138
0	16	10	7	135	170	153
0	16	11	8	135	170	153
0	16	13	10	135	170	153
0	16	15	12	135	170	153
0/3/6	16	10	7	160	195	178
0/3/6	16	11	8	160	195	178
0/3/6	16	13	10	160	195	178
0/3/6	16	15	12	160	195	178
0/3/6	16	11	8	180	215	198
0/3/6	16	13	10	180	215	198
0/3/6	16	15	12	180	215	198
0	16	11	8	200	235	218
0	16	13	10	200	235	218
0	16	15	12	200	235	218
0/3/6	16	11	8	240	275	258
0/3/6	16	13	10	240	275	258
0/3/6	16	15	12	240	275	258
0	16	11	8	280	315	298
0	16	13	10	280	315	298
0	16	15	12	280	315	298

LinkSymphoKnee – Centralizers




Ø mm	Length (L) mm
12	15
14	15
16	15
18	15
20	15
22	15
24	15

LinkSymphoKnee Femoral Augments – Distal




cemented		
Size	Side	Height mm
0	Medial-Right/Lateral-Left	5
0	Lateral-Right/Medial-Left	5
1-2	Medial-Right/Lateral-Left	5
1-2	Lateral-Right/Medial-Left	5
3-4	Medial-Right/Lateral-Left	5
3-4	Medial-Right/Lateral-Left	10
3-4	Lateral-Right/Medial-Left	5
3-4	Lateral-Right/Medial-Left	10
5-6	Medial-Right/Lateral-Left	5
5-6	Medial-Right/Lateral-Left	10
5-6	Lateral-Right/Medial-Left	5
5-6	Lateral-Right/Medial-Left	10
7-8	Medial-Right/Lateral-Left	5
7-8	Medial-Right/Lateral-Left	10
7-8	Lateral-Right/Medial-Left	5
7-8	Lateral-Right/Medial-Left	10
9-10	Medial-Right/Lateral-Left	5
9-10	Medial-Right/Lateral-Left	10
9-10	Lateral-Right/Medial-Left	5
9-10	Lateral-Right/Medial-Left	10

LinkSymphoKnee Femoral Augments – Posterior



cemented		
Size	Side	Height mm
0	Medial-Right/Lateral-Left	5
0	Lateral-Right/Medial-Left	5
1-2	Medial-Right/Lateral-Left	5
1-2	Lateral-Right/Medial-Left	5
3-4	Medial-Right/Lateral-Left	5
3-4	Medial-Right/Lateral-Left	10
3-4	Lateral-Right/Medial-Left	5
3-4	Lateral-Right/Medial-Left	10
5-6	Medial-Right/Lateral-Left	5
5-6	Medial-Right/Lateral-Left	10
5-6	Lateral-Right/Medial-Left	5
5-6	Lateral-Right/Medial-Left	10
7-8	Medial-Right/Lateral-Left	5
7-8	Medial-Right/Lateral-Left	10
7-8	Lateral-Right/Medial-Left	5
7-8	Lateral-Right/Medial-Left	10
9-10	Medial-Right/Lateral-Left	5
9-10	Medial-Right/Lateral-Left	10
9-10	Lateral-Right/Medial-Left	5
9-10	Lateral-Right/Medial-Left	10


LinkSymphoKnee Femoral Augments – L-Shaped



cemented

Size	Side	Height mm
0	Medial-Right/Lateral-Left	10
0	Lateral-Right/Medial-Left	10
1-2	Medial-Right/Lateral-Left	10
1-2	Lateral-Right/Medial-Left	10
3-4	Medial-Right/Lateral-Left	15
3-4	Lateral-Right/Medial-Left	15
5-6	Medial-Right/Lateral-Left	15
5-6	Lateral-Right/Medial-Left	15
7-8	Medial-Right/Lateral-Left	15
7-8	Lateral-Right/Medial-Left	15
9-10	Medial-Right/Lateral-Left	15
9-10	Lateral-Right/Medial-Left	15

LinkSymphoKnee Tibial Augments



cemented

Size	Side	Height mm
1-2	Medial-Right/Lateral-Left	5
1-2	Medial-Right/Lateral-Left	10
1-2	Medial-Right/Lateral-Left	15
1-2	Lateral-Right/Medial-Left	5
1-2	Lateral-Right/Medial-Left	10
1-2	Lateral-Right/Medial-Left	15
3-4	Medial-Right/Lateral-Left	5
3-4	Medial-Right/Lateral-Left	10
3-4	Medial-Right/Lateral-Left	15
3-4	Lateral-Right/Medial-Left	5
3-4	Lateral-Right/Medial-Left	10
3-4	Lateral-Right/Medial-Left	15
5-6	Medial-Right/Lateral-Left	5
5-6	Medial-Right/Lateral-Left	10
5-6	Medial-Right/Lateral-Left	15
5-6	Lateral-Right/Medial-Left	5
5-6	Lateral-Right/Medial-Left	10
5-6	Lateral-Right/Medial-Left	15
7-8	Medial-Right/Lateral-Left	5
7-8	Medial-Right/Lateral-Left	10
7-8	Medial-Right/Lateral-Left	15
7-8	Lateral-Right/Medial-Left	5
7-8	Lateral-Right/Medial-Left	10
7-8	Lateral-Right/Medial-Left	15
9-10	Medial-Right/Lateral-Left	5
9-10	Medial-Right/Lateral-Left	10
9-10	Medial-Right/Lateral-Left	15
9-10	Lateral-Right/Medial-Left	5
9-10	Lateral-Right/Medial-Left	10
9-10	Lateral-Right/Medial-Left	15

- 1 Koninckx, A., Deltour, A. & Thienpont, E. Femoral sizing in total knee arthroplasty is rotation dependant. *Knee Surg Sports Traumatol Arthrosc* 22, 2941–2946 (2014). <https://doi.org/10.1007/s00167-013-2707-5>
- 2 Kia M, Wright TM, Cross MB, Mayman DJ, Pearle AD, Sculco PK, Westrich GH, Imhauser CW. Femoral Component External Rotation Affects Knee Biomechanics: A Computational Model of Posterior-stabilized TKA. *Clin Orthop Relat Res*. 2018 Jan
- 3 Chadd W. Clary, Clare K. Fitzpatrick, Lorin P. Maletsky, Paul J. Rullkoetter, The influence of total knee arthroplasty geometry on mid-flexion stability: An experimental and finite element study, *Journal of Biomechanics*, Volume 46, Issue 7, 2013, Pages 1351-1357, ISSN 0021-9290, <https://doi.org/10.1016/j.jbiomech.2013.01.025>. (<https://www.sciencedirect.com/science/article/pii/S0021929013000675>)
- 4 Graichen H. TKA revision - reasons, challenges and solutions. *J Orthop*. 2014 Mar 26;11(1):1-4. doi: 10.1016/j.jor.2014.01.005. PMID: 24719528; PMCID: PMC3978737
- 5 Petersen W, Rembitzki IV, Brüggemann GP, Ellermann A, Best R, Koppenburg AG, Liebau C. Anterior knee pain after total knee arthroplasty: a narrative review. *Int Orthop*. 2014 Feb;38(2):319-28. doi: 10.1007/s00264-013-2081-4. Epub 2013 Sep 22. PMID: 24057656; PMCID: PMC3923935
- 6 Meijerink HJ, Barink M, van Loon CJ, Schwering PJ, Donk RD, Verdonschot N, de Waal Malefijt MC. The trochlea is medialized by total knee arthroplasty: an intraoperative assessment in 61 patients. *Acta Orthop*. 2007 Feb;78(1):123-7. doi: 10.1080/17453670610013529. PMID: 17453403.
- 7 Completo A, Simões JA, Fonseca F, Oliveira M. The influence of different tibial stem designs in load sharing and stability at the cement-bone interface in revision TKA. *Knee*. 2008 Jun;15(3):227-32
- 8 Dennis DA, Komistek RD. Mobilebearing total knee arthroplasty: design factors in minimizing wear. *Clin Orthop Relat Res*. 2006 Nov;452:70-7.
- 9 Walker PS, Hsu HP, Zimmerman RA. A comparative study of uncemented tibial components. *J Arthroplasty*. 1990 Sep;5(3):245-53.
- 10 Stem Fixation in Revision Total Knee Arthroplasty: Indications, Stem Dimensions, and Fixation Methods - PMC (nih.gov)
- 11 The anatomy of the tibial intramedullary canal - PubMed (nih.gov)
- 12 Barrack, R.L.; Rorabeck, C.; Burt, M.; Sawhney, J., 1999; Barrack, R.L.; Stanley, T.; Burt, M.; Hopkins, S., 2004.
- 13 Ries MD, Suzuki Y, Renowitzky G, Lotz JC, Barrack RL, Bourne RB, Rorabeck CH. Effect of cementless bowed stem distal surface contour and coronal slot on femoral bone strains and torsional stability. *J Arthroplasty*. 2003 Jun;18(4):494-8. doi: 10.1016/s0883-5403(03)00069-x. PMID: 12820094.
- 14 Polizzotti G, Lamberti A, Mancino F, Baldini A. New Horizons of Cementless Total Knee Arthroplasty. *J Clin Med*. 2023 Dec 30;13(1):233. doi: 10.3390/jcm13010233. PMID: 38202240; PMCID: PMC10780266.
- 15 Gross TP, Liu F. Total knee arthroplasty with fully porous-coated stems for the treatment of large bone defects. *J Arthroplasty*. 2013 Apr;28(4):598-603. doi: 10.1016/j.arth.2012.07.021. Epub 2012 Nov 8. PMID: 23141862.
- 16 Thienpont E. Titanium niobium nitride knee implants are not inferior to chrome cobalt components for primary total knee arthroplasty at medium-term follow-up. *Arch Orthop Trauma Surg*. 2023 Aug;143(8):5269-5275. doi: 10.1007/s00402-022-04754-1. Epub 2023 Jan 3. PMID: 36595031.
- 17 Bader R, Bergschmidt P, Fritsche A, Ansorge S, Thomas P, Mittelmeier W. Alternative Werkstoffe und Lösungen in der Knieendoprothetik für Patienten mit Metallallergie. *Orthopäde*, 2008; 37:136-142
- 18 S. M. Kurtz, „The Origins and Adaptations of UHMWPE for Knee replacement“, in *UHMWPE Biomaterials Handbook*, S. M. Kurtz, d., Burlington, MA Academic Press 2009.
- 19 S. M. Kurtz, „Advances in th the processing, sterilization, and crosslinking of ultra-high molecular weight polyethylene for total joint arthroplasty“, *Biomaterials* 1999; 20:1659-1687.
- 20 E. M. Brach del Prever, „UHMWPE for arthroplasty: past ot future?“, *J Orthopaed Traumatol* 2009;10:1-8.
- 21 Bracco, P., & Oral, E. (2011). Vitamin E-stabilized UHMWPE for total joint implants: a review. *Clinical orthopaedics and related research*, 469(8), 2286–2293. <https://doi.org/10.1007/s11999-010-1717-6>
- 22 S. M. Kurtz, „The Clinical Performance of UHMWPE in Knee Replacement“, in *UHMWPE Biomaterials Handbook*, S. M. Kurtz, Ed., Burlington, MA Academic Press 2009.
- 23 M. C. Forster, „Survival Analysis of Primary Cemented Total Knee Arthroplasty“, *J Arthroplasty* 2003; 18:265-270.
- 24 Swedish Knee Arthroplasty Register, Annual Report 2015, www.myknee.se
- 25 National Joint Registry for England, Wales, Northern Ireland and the Isle of Man, Annual Report 2015, www.njrreports.org.uk
- 26 Oral, E., Christensen, S. D., Malhi, A. S., Wannomae, K. K., & Muratoglu, O. K. (2006). Wear resistance and mechanical properties of highly cross-linked, ultrahigh-molecular weight polyethylene doped with vitamin E. *The Journal of arthroplasty*, 21(4), 580–591. <https://doi.org/10.1016/j.arth.2005.07.009>
- 27 E. Oral, „Characterization of Irradiated Blends of Alpha-tocopherol and UHMWPE“, *Biomaterials* 2005; 26(33):6657-6663.
- 28 E. Oral, „Alpha-tocopherol-doped Irradiated UHMWPE for High Fatigue Resistance and Low Wear“, *Biomaterials* 2004; 25:5515-5522
- 29 O. K. Muratoglu, „Knee Simulator Testing of Conventional and Cross-Linked Polyethylene Tibial Inserts“, *J Arthroplasty* 2004; 19:887-897.
- 30 O. K. Muratoglu, „Aggressive Wear Testing of a Cross-Linked Polyethylene in Total Knee Arthroplasty“, *Clin Orthop Relat Res*. 2002; 404:89-95.
- 31 J. T. Hodrick, „Highly Crosslinked Polyethylene is Safe for Use in Total Knee Arthroplasty“, *Clin Orthop Relat Res*. 2008; 466:2806-2812.
- 32 R. N. De Steiger, „Lower Prosthesis-specific 10-year Revision Rate with Crosslinked than with Non-crosslinked Polyethylene in Primary Total Knee Arthroplasty“, *Acta Orthopædica* 2015; 86 (6):721-727.
- 33 E. Oral, „Highly Crosslinked UHMWPE Doped with Vitamin E“, in *UHMWPE Biomaterials Handbook*, S. M. Kurtz, Ed., Burlington, MA Academic Press 2009.
- 34 S. M. Kurtz, „Vitamin-E-Blended UHMWPE Biomaterials“, in *UHMWPE Biomaterials Handbook*, S. M. Kurtz, Ed., Burlington, MA Academic Press 2009.
- 35 S. M. Kurtz, „Trace Concentration of Vitamin E Protect Radiation Crosslinked UHMWPE from Oxidative Degradation“, *J Biomed Mater Res A* 2008; 549-563
- 36 Gigante, Antonio, et al. "Effectiveness of Vitamin-E-Doped Polyethylene in Joint Replacement: A Literature Review." *PubMed Central (PMC)*, 8 Sept. 2015, www.ncbi.nlm.nih.gov/pmc/articles/PMC4598683

Waldemar Link GmbH & Co. KG, Hamburg

All content in this catalog, including text, pictures and data, is protected by law. Every instance of use, whether in part or in whole and which is not permitted by law, is subject to our prior consent. In particular, this applies to the reproduction, editing, translation, publishing, saving, processing, or passing on of content stored in databases or other electronic media and systems, in any manner or form. The information in the catalogs is solely intended to describe the products and does not constitute a guarantee.

The Surgical Technique described has been written to the best of our knowledge and belief, but it does not relieve the surgeon of his/her responsibility to duly consider the particularities of each individual case.

Products shown in this document may not be available in your country. The product availability is subject to the approval and/or registration regulations of the respective country. Please contact Waldemar Link GmbH & Co. KG if you have questions about the availability of LINK products in your country.

Waldemar Link GmbH & Co. KG and/or other corporate affiliated entities own, use or have applied for the following trademarks in many jurisdictions: LINK, BiMobile, SP II, Modell Lubinus, E-Dur, EndoDur, T.O.P. II, BetaCup, CombiCup PF, CombiCup SC, CombiCup R, MobileLink, C.F.P., LCU, SP-CL, LCP, MIT-H, Endo-Model, Endo-Model SL, MP, MEGASYSTEM-C, GEMINI SL, SPAR-K, LCK, Link OptiStem, HX, TiCaP, X-LINKed, PorAg, LINK PorEx, BiPorEx, PorEx-Z, TrabecuLink, Tilastan, customLINK, RescueSleeve, Stactip, VACUCAST.

Other trademarks and trade names may be used in this document to refer to either the entities claiming the marks and/or names or their products and are the property of their respective owners.



For more information please register for our LINK Media Library (link-ortho.com)

 Waldemar Link GmbH & Co. KG

Barkhausenweg 10 · 22339 Hamburg · Germany
Phone +49 40 53995-0 · info@link-ortho.com
www.link-ortho.com

